

## Table of contents

\\Study Protocols	
BRAIN	
Stroke	
Investigate - E161686 v1b	
AAHead_Scout_32ch-head-coil	
CVR_pace_2.5mm_iso_v1_TEST_SCAN	
CVR_pace_2.5mm_iso_v1	
AAHead_Scout_32ch-head-coil	
t1_mprage_sag_p3_iso_Munich	
t2_3D_FLAIR	v2
t2_space_v4	
PD_VIBE_v1a	
SWI_v6c	
t2_blade_dark-fluid_tra_v2b	
DTI_custom_v1_AP	
DTI_custom_v1_rev	
localizer	
TOF	
PC_Carotids	
localizer	(C2-C3)
PC_SACSF	
localizer	
vessel_scout	
PC_Sinus	
t1_fl2d_sag_localiser	
PC_Aqueduct	
AAHead_Scout_32ch-head-coil	
IR-SPGR_TI_600_v4	
IR-SPGR_TI_1500_v4	
IR-SPGR_IROff_2_v4	
IR-SPGR_IROff_5_v4	
IR-SPGR_IROff_12_v4	
t1_vibe_sag_DCE_2mm	

\\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\AAHead\_Scout\_32ch-head-coil

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3.15 ms
TE	1.37 ms
Flip angle	8 deg

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

**Resolution - iPAT**

Reference scan mode	Integrated
---------------------	------------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off

**System - Miscellaneous**

Coil Select Mode	Default
------------------	---------

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Flip angle	8 deg
Measurements	1
Time to center	6.2 s

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1

**Inline - MapIt**

Contrasts	1
TR	3.15 ms
TE	1.37 ms

**Sequence - Part 1**

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
------	-----

\\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\CVR\_pace\_2.5mm\_iso\_v1\_TEST\_SCAN

TA: 0:12 PM: REF Voxel size: 2.5×2.5×2.5 mmPAT: 2 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	On
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	52
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	235 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
TR	3000 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3000 ms
TE	30.0 ms
MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	235 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
Base resolution	94
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

Accel. mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24

**Resolution - iPAT**

Reference scan mode	EPI/separate
---------------------	--------------

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	On

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	52
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	235 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
TR	3000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain

**System - Miscellaneous**

Coil Select Mode	On - AutoCoilSelect
------------------	---------------------

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Performance

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	235 mm
R >> L	235 mm
F >> H	130 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3000 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion correction	Off
Spatial filter	Off
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.52 ms
Bandwidth	2312 Hz/Px

**Sequence - Part 2**

EPI factor	94
------------	----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\CVR\_pace\_2.5mm\_iso\_v1

TA: 12:30 PM: REF Voxel size: 2.5×2.5×2.5 mmPAT: 2 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	On
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	52
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
FoV read	235 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
TR	3000 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3000 ms
TE	30.0 ms
MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	247
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	235 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
Base resolution	94
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

Accel. mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24

**Resolution - iPAT**

Reference scan mode	EPI/separate
---------------------	--------------

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	On

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	52
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	235 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
TR	3000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	Head > Brain

**System - Miscellaneous**

Coil Select Mode	On - AutoCoilSelect
------------------	---------------------

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Performance

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	235 mm
R >> L	235 mm
F >> H	130 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3000 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	3
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Active
Motion correction	Off
Spatial filter	Off
Measurements	247
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.52 ms
Bandwidth	2312 Hz/Px

**Sequence - Part 2**

EPI factor	94
------------	----

\\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\AAHead\_Scout\_32ch-head-coil

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3.15 ms
TE	1.37 ms
Flip angle	8 deg

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

**Resolution - iPAT**

Reference scan mode	Integrated
---------------------	------------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off



**System - Miscellaneous**

Coil Select Mode	Default
------------------	---------

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Flip angle	8 deg
Measurements	1
Time to center	6.2 s

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1

**Inline - MapIt**

Contrasts	1
TR	3.15 ms
TE	1.37 ms

**Sequence - Part 1**

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\t1\_mprage\_sag\_p3\_iso\_Munich

TA: 3:45 PM: FIX Voxel size: 1.0×1.0×1.0 mmPAT: 3 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2500.0 ms
TE	4.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2500.0 ms
TE	4.37 ms
Magn. preparation	Non-sel. IR
T1	1100 ms
Flip angle	7 deg
Fat suppr.	Water excit. fast
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	7/8
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	2500.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	256 mm
F >> H	256 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2500.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI	1100 ms
Fat suppr.	Water excit. fast
Dark blood	Off
FoV read	256 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	7 deg
Measurements	1
TR	2500.0 ms
TE	4.37 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	11.1 ms
Bandwidth	140 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	192

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b1t2\_3D\_FLAIR v2

TA: 5:57 PM: REF Voxel size: 1.0×1.0×1.0 mmPAT: 3 Rel. SNR: 1.00 : spcir

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	5000 ms
TE	388 ms
Averages	1.0
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	5000 ms
TE	388 ms
MTC	Off
Magn. preparation	Non-sel. T2-IR
TI 1	1800 ms
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Blood suppr.	Off
Restore magn.	Off

**Contrast - Dynamic**

Averages	1.0
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Allowed
Slice partial Fourier	Off

**Resolution - Common**

Interpolation	Off
---------------	-----

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	5000 ms
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Restore magn.	Off
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
-------------------	-----

**Geometry - Tim Planning Suite**

Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Sagittal
Rotation	0.00 deg
A >> P	256 mm
F >> H	256 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
Trigger delay	0 ms
TR	5000 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. T2-IR
TI 1	1800 ms
Fat suppr.	Fat sat.
Dark blood	Off
FoV read	256 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Flow comp.	No
Echo spacing	3.46 ms
Adiabatic-mode	Off
Bandwidth	751 Hz/Px

**Sequence - Part 2**

Echo train duration	865 ms
RF pulse type	Normal
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
Turbo factor	278

**Sequence - Assistant**

Allowed delay	30 s
---------------	------

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\t2\_space\_v4

TA: 3:42 PM: REF Voxel size: 0.9×0.9×0.9 mmPAT: 4 Rel. SNR: 1.00 : spcR

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	18.2 %
Slices per slab	176
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	0.90 mm
TR	3200 ms
TE	408 ms
Averages	1.4
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3200 ms
TE	408 ms
MTC	Off
Magn. preparation	None
Fat suppr.	None
Blood suppr.	Off
Restore magn.	On

**Contrast - Dynamic**

Averages	1.4
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	240 mm
FoV phase	100.0 %
Slice thickness	0.90 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Allowed
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	2
Ref. lines 3D	24
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	18.2 %
Slices per slab	176
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	0.90 mm
TR	3200 ms
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	90.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Restore magn.	On
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Performance
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
Trigger delay	0 ms
TR	3200 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	240 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
----------	-----

**Inline - Common**

Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Flow comp.	No
Echo spacing	3.61 ms
Adiabatic-mode	Off
Bandwidth	723 Hz/Px

**Sequence - Part 2**

Echo train duration	910 ms
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
Flip angle mode	T2 var
Turbo factor	282

**Sequence - Assistant**

Allowed delay	30 s
---------------	------

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\PD\_VIBE\_v1a

TA: 1:57 PM: FIX Voxel size: 1.0×1.0×1.0 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A5.1 F8.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	6.04 ms
TE	2.44 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP;SP1

**Contrast - Common**

TR	6.04 ms
TE	2.44 ms
Flip angle	2.0 deg
Fat suppr.	None
Water suppr.	None
Dixon	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
View sharing	Off

**Resolution - Common**

Interpolation	Off
---------------	-----

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A5.1 F8.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	192
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	6.04 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	L0.0 A5.1 F8.3 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	L0.0 A5.1 F8.3
R	0.0 mm
A	5.1 mm
F	8.3 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	None
Water suppr.	None
Dixon	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
-------------------	-----



**Geometry - Tim Planning Suite**

Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	Head > Basis
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.0 A5.1 F8.3 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	256 mm
F >> H	256 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	2.0 deg
Measurements	1
Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
Time to center	58.9 s

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	2.0 deg
Measurements	1
Contrasts	1
TR	6.04 ms
TE	2.44 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	1
Optimization	Min. TE
Multi-slice mode	Sequential
Bandwidth	200 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\SWI\_v6c

TA: 4:02 PM: FIX Voxel size: 0.6×0.6×3.0 mmPAT: 2 Rel. SNR: 1.00 : swi\_r

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A2.3 H2.2 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	23.1 %
Slices per slab	52
FoV read	240 mm
FoV phase	81.3 %
Slice thickness	3.00 mm
TR	28.0 ms
TE	20.00 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	28.0 ms
TE	20.00 ms
MTC	Off
Magn. preparation	None
Flip angle	9 deg
Fat suppr.	None
Water suppr.	None
SWI	On

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	240 mm
FoV phase	81.3 %
Slice thickness	3.00 mm
Base resolution	384
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off

**Resolution - Common**

Interpolation	Off
---------------	-----

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A2.3 H2.2 mm
Orientation	Transversal
Phase enc. dir.	R >> L
Slice oversampling	23.1 %
Slices per slab	52
FoV read	240 mm
FoV phase	81.3 %
Slice thickness	3.00 mm
TR	28.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	L0.0 A2.3 H2.2 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	L0.0 A2.3 H2.2
L	0.0 mm
A	2.3 mm
H	2.2 mm
Initial Rotation	89.61 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H

**Geometry - Tim Planning Suite**

Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.0 A2.3 H2.2 mm
Orientation	Transversal
Rotation	89.61 deg
R >> L	195 mm
A >> P	240 mm
F >> H	156 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slab-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	28.0 ms
Concatenations	1
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	240 mm
FoV phase	81.3 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	9 deg
Measurements	1
Contrasts	1
TR	28.0 ms
TE	20.00 ms

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	On
Phase stabilisation	Off
Asymmetric echo	Off
Contrasts	1
Flow comp.	Yes
Multi-slice mode	Interleaved
Bandwidth	120 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	None
RF pulse type	Fast
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	30 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b1t2\_blade\_dark-fluid\_tra\_v2b

TA: 3:22 PM: REF Voxel size: 0.9×0.9×3.0 mmPAT: 2 Rel. SNR: 1.00 : tirB\_rr

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	52
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0.0 %
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10000.0 ms
TE	124 ms
Averages	1
Concatenations	2
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	10000.0 ms
TE	124 ms
TD	0.0 ms
MTC	Off
Magn. preparation	Slice-sel. IR
T1	2606 ms
Flip angle	130 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Water suppr.	None
Restore magn.	Off
Freeze suppressed tissue	On

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
Base resolution	256
BLADE coverage	100.0 %
Trajectory	BLADE

**Resolution - Common**

Interpolation	Off
---------------	-----

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	8
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	52
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	240 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	2

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	90.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Water suppr.	None
Restore magn.	Off
Special sat.	Parallel F
Gap	10 mm
Thickness	70 mm

**Geometry - Navigator**

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
R >> L	240 mm
A >> P	240 mm
F >> H	156 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	10000.0 ms
Concatenations	2

**Physio - Cardiac**

Magn. preparation	Slice-sel. IR
TI	2606 ms
Fat suppr.	Fat sat.
Dark blood	Off
FoV read	240 mm
FoV phase	100.0 %
BLADE coverage	100.0 %
Trajectory	BLADE

**Physio - PACE**

Resp. control	Off
Concatenations	2

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Compensate T2 decay	Off
Contrasts	1
Flow comp.	Read
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	8.26 ms
Bandwidth	362 Hz/Px

**Sequence - Part 2**

Define	Turbo factor
Echo trains per slice	9
Phase correction	Automatic
Acoustic noise reduction	None
RF pulse type	Low SAR
Gradient mode	Fast
Hyperecho	Off
WARP	Off
Motion correction	On
Red. EC sensitivity	Off
Turbo factor	28

**Sequence - Assistant**

Mode	Min flip angle
Min flip angle	130 deg
Allowed delay	30 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\DTI\_custom\_v1\_AP

TA: 11:12 PM: REF Voxel size: 2.0×2.0×2.0 mmPAT: 4 Rel. SNR: 1.00 : epse

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4300 ms
TE	74.0 ms
Averages	1
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	4300 ms
TE	74.0 ms
MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Strong

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
-------------	--------------

**Resolution - iPAT**

Accel. factor PE	2
Ref. lines PE	40
Accel. factor slice	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	On
Dynamic Field Corr.	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4300 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.7 P3.0 H31.9
L	0.7 mm
P	3.0 mm
H	31.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Matrix Optimization	Performance
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Rotation	0.00 deg
A >> P	256 mm
R >> L	256 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Standard

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	4300 ms
Concatenations	1

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Diff - Neuro**

Diffusion mode	Free
Diff. directions	151
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	2000 s/mm <sup>2</sup>
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	40

**Diff - Body**

Diffusion mode	Free
Diff. directions	151
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	2000 s/mm <sup>2</sup>
b-value	1
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm <sup>2</sup>
Noise level	40

**Diff - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Optimization	None
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.78 ms
Bandwidth	1446 Hz/Px

**Sequence - Part 2**

EPI factor	128
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard

**Sequence - pTX Pulses**

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\DTI\_custom\_v1\_rev

TA: 0:35 PM: REF Voxel size: 2.0×2.0×2.0 mmPAT: 4 Rel. SNR: 1.00 : epse

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Phase enc. dir.	P >> A
AutoAlign	---
Phase oversampling	0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4300 ms
TE	74.0 ms
Concatenations	1
Filter	Raw filter, Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	4300 ms
TE	74.0 ms
MTC	Off
Magn. preparation	None
Fat suppr.	Fat sat.
Fat sat. mode	Strong

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	2
Ref. lines PE	40

**Resolution - iPAT**

Accel. factor slice	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	On
Dynamic Field Corr.	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	72
Dist. factor	0 %
Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Phase enc. dir.	P >> A
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	4300 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Phase enc. dir.	P >> A
AutoAlign	---
Initial Position	L0.7 P3.0 H31.9
L	0.7 mm
P	3.0 mm
H	31.9 mm
Initial Rotation	180.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Fat sat. mode	Strong
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L



**System - Miscellaneous**

Coronal	A >> P
Transversal	H >> F
Coil Combine Mode	Adaptive Combine
Matrix Optimization	Performance
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.7 P3.0 H31.9 mm
Orientation	Transversal
Rotation	180.00 deg
A >> P	256 mm
R >> L	256 mm
F >> H	144 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Standard

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	4300 ms
Concatenations	1

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Diff - Neuro**

Diffusion mode	MDDW
Diff. directions	6
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	0 s/mm <sup>2</sup>
b-value	3
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
FA maps	Off
Mosaic	Off
Tensor	Off
Noise level	40

**Diff - Body**

Diffusion mode	MDDW
----------------	------

**Diff - Body**

Diff. directions	6
Diffusion Scheme	Monopolar
Diff. weightings	1
b-value	0 s/mm <sup>2</sup>
b-value	3
Diff. weighted images	On
Trace weighted images	Off
ADC maps	Off
Exponential ADC Maps	Off
FA maps	Off
Invert Gray Scale	Off
Calculated Image	Off
b-Value >=	0 s/mm <sup>2</sup>
Noise level	40

**Diff - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Optimization	None
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.78 ms
Bandwidth	1446 Hz/Px

**Sequence - Part 2**

EPI factor	128
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard

**Sequence - pTX Pulses**

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\localizer

TA: 0:11 PM: ISO Voxel size: 1.2×1.2×10.0 mmPAT: Off Rel. SNR: 1.00 : fl\_rr

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	2
Dist. factor	70 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	9.7 ms
TE	3.80 ms
Averages	1
Concatenations	8
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	9.7 ms
TE	3.80 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	40 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Contrast - Dynamic**

Multiple series	Each measurement
-----------------	------------------

**Resolution - Common**

FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
Base resolution	256
Phase resolution	50 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	2
Dist. factor	70 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	9.7 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	8

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P

**Geometry - AutoAlign**

Slice group	2
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	ISO
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	9.7 ms
Concatenations	8
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	300 mm
FoV phase	100.0 %
Phase resolution	50 %

**Physio - PACE**

Resp. control	Off
Concatenations	8

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	40 deg
Measurements	1
Contrasts	1
TR	9.7 ms
TE	3.80 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	Read
Multi-slice mode	Sequential
Bandwidth	320 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\TOF

TA: 2:45 PM: ISO Voxel size: 0.5×0.5×0.8 mmPAT: 2 Rel. SNR: 1.00 : fl\_r

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	6
Dist. factor	-20.83 %
Position	R4.6 A0.9 F54.5 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	16.7 %
Slices per slab	24
FoV read	200 mm
FoV phase	75.5 %
Slice thickness	0.80 mm
TR	20.0 ms
TE	3.51 ms
Averages	1
Concatenations	6
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	HE3,4;NE1,2;SP1

**Contrast - Common**

TR	20.0 ms
TE	3.51 ms
TD	0.000 ms
MTC	Off
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	200 mm
FoV phase	75.5 %
Slice thickness	0.80 mm
Base resolution	384
Phase resolution	70 %
Slice resolution	50 %
Phase partial Fourier	7/8
Slice partial Fourier	7/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	6
Dist. factor	-20.83 %
Position	R4.6 A0.9 F54.5 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice oversampling	16.7 %
Slices per slab	24
FoV read	200 mm
FoV phase	75.5 %
Slice thickness	0.80 mm
TR	20.0 ms
Multi-slice mode	Sequential
Series	Descending
Concatenations	6

**Geometry - AutoAlign**

Slab group	1
Position	R4.6 A0.9 F54.5 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.6 A0.9 F54.5
R	4.6 mm
A	0.9 mm
F	54.5 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	None
Water suppr.	None
Special sat.	Tracking H
Gap	10 mm
Thickness	40 mm

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
-------------------	-----

**Geometry - Tim Planning Suite**

Table position	F
Table position	55 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	ISO
Table position	F
Table position	55 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	20.0 ms
Concatenations	6

**Physio - Cardiac**

Fat suppr.	None
Dark blood	Off
FoV read	200 mm
FoV phase	75.5 %
Phase resolution	70 %

**Angio - Common**

TONE ramp	80 %
Flow direction	F >> H

**Angio - Common**

Flip angle	20 deg
MTC	Off
Measurements	1
3D centric reordering	On

**Angio - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Angio - MIP**

MIP-Sag	On
MIP-Cor	On
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Angio - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	Yes
Multi-slice mode	Sequential
Bandwidth	250 Hz/Px

**Sequence - Part 2**

Gradient mode	Fast
RF spoiling	On

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\PC\_Carotids

TA: 1:39 PM: ISO Voxel size: 1.0×1.0×5.0 mmPAT: 2 Rel. SNR: 1.00 : fl\_r

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	R4.9 A18.6 F16.5 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	19.60 ms
TE	5.82 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	HE3,4;NE1,2

**Contrast - Common**

TR	19.60 ms
TE	5.82 ms
TD	0 ms
Flip angle	12 deg
Wrap-up Magn.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	160
Phase resolution	100 %
Phase partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
----------	--------

**Resolution - iPAT**

Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
POCS	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	R4.9 A18.6 F16.5 mm
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	19.60 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.9 A18.6 F16.5 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	R4.9 A18.6 F16.5
R	4.9 mm
A	18.6 mm
F	16.5 mm
Initial Rotation	90.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Wrap-up Magn.	None
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	F
Table position	17 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	ISO
Table position	F
Table position	17 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.9 A18.6 F16.5 mm
Orientation	Transversal
Rotation	90.00 deg
R >> L	160 mm
A >> P	160 mm
F >> H	5 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	Pulse/Retro
Average cycle	No Signal ms
Average cycle	No Signal ms
Calculated phases	32
TR	19.60 ms
Concatenations	1
Segments	1
Arrhythmia detection	None

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Angio - Common**

Flow mode	Single dir.
Encodings	1
Velocity enc.	70 cm/s
Direction	Through plane
Rephased images	On
Magnitude images	On

**Angio - Common**

Magnitude sum	Off
Phase images	On

**Angio - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Angio - Cardiac**

Inline Evaluation	Off
TE	5.82 ms
TR	19.60 ms
Save original images	On

**Angio - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Angio - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Flow comp.	Yes
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	9.8 ms
Sequence type	Gre
Bandwidth	202 Hz/Px

**Sequence - Part 2**

Define	Segments
Segments	1
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s



## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\localizer (C2-C3)

TA: 0:11 PM: ISO Voxel size: 1.2×1.2×10.0 mmPAT: Off Rel. SNR: 1.00 : fl\_rr

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	2
Dist. factor	70 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	9.7 ms
TE	3.80 ms
Averages	1
Concatenations	8
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	HE1-4;NE1,2

**Contrast - Common**

TR	9.7 ms
TE	3.80 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	40 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Contrast - Dynamic**

Multiple series	Each measurement
-----------------	------------------

**Resolution - Common**

FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
Base resolution	256
Phase resolution	50 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	2
Dist. factor	70 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	9.7 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	8

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P

**Geometry - AutoAlign**

Slice group	2
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	ISO
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	9.7 ms
Concatenations	8
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	300 mm
FoV phase	100.0 %
Phase resolution	50 %

**Physio - PACE**

Resp. control	Off
Concatenations	8

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	40 deg
Measurements	1
Contrasts	1
TR	9.7 ms
TE	3.80 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	Read
Multi-slice mode	Sequential
Bandwidth	320 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\PC\_SACSF

TA: 1:55 PM: ISO Voxel size: 0.8×0.8×5.0 mmPAT: 2 Rel. SNR: 1.00 : fl\_r

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	L0.1 A7.3 F11.2 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	25.18 ms
TE	8.45 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	HE3,4;NE1,2

**Contrast - Common**

TR	25.18 ms
TE	8.45 ms
TD	0 ms
Flip angle	12 deg
Wrap-up Magn.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	192
Phase resolution	100 %
Phase partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On
POCS	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	L0.1 A7.3 F11.2 mm
Orientation	Transversal
Phase enc. dir.	R >> L
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	25.18 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.1 A7.3 F11.2 mm
Orientation	Transversal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L0.1 A7.3 F11.2
L	0.1 mm
A	7.3 mm
F	11.2 mm
Initial Rotation	90.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Wrap-up Magn.	None
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	F
Table position	11 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	ISO
Table position	F
Table position	11 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.1 A7.3 F11.2 mm
Orientation	Transversal
Rotation	90.00 deg
R >> L	160 mm
A >> P	160 mm
F >> H	5 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	Pulse/Retro
Average cycle	No Signal ms
Average cycle	No Signal ms
Calculated phases	32
TR	25.18 ms
Concatenations	1
Segments	1
Arrhythmia detection	None

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Angio - Common**

Flow mode	Single dir.
Encodings	1
Velocity enc.	6 cm/s
Direction	Through plane
Rephased images	On
Magnitude images	On

**Angio - Common**

Magnitude sum	Off
Phase images	On

**Angio - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Angio - Cardiac**

Inline Evaluation	Off
TE	8.45 ms
TR	25.18 ms
Save original images	On

**Angio - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Angio - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Flow comp.	Yes
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	12.6 ms
Sequence type	Gre
Bandwidth	200 Hz/Px

**Sequence - Part 2**

Define	Segments
Segments	1
RF pulse type	Normal
Gradient mode	Fast*
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\localizer

TA: 0:11 PM: ISO Voxel size: 1.2×1.2×10.0 mmPAT: Off Rel. SNR: 1.00 : fl\_rr

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	2
Dist. factor	70 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	9.7 ms
TE	3.80 ms
Averages	1
Concatenations	8
Filter	Distortion Corr.(2D), Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	9.7 ms
TE	3.80 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	40 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Contrast - Dynamic**

Multiple series	Each measurement
-----------------	------------------

**Resolution - Common**

FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
Base resolution	256
Phase resolution	50 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	300 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	2
Dist. factor	70 %
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	9.7 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	8

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P

**Geometry - AutoAlign**

Slice group	2
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Position	L0.0 A30.0 H0.0 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	ISO
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	9.7 ms
Concatenations	8
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	300 mm
FoV phase	100.0 %
Phase resolution	50 %

**Physio - PACE**

Resp. control	Off
Concatenations	8

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	40 deg
Measurements	1
Contrasts	1
TR	9.7 ms
TE	3.80 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	Read
Multi-slice mode	Sequential
Bandwidth	320 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s



**\\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\vessel\_scout**

TA: 0:25 PM: FIX Voxel size: 1.4×1.4×5.0 mmPAT: 2 Rel. SNR: 1.00 : pc

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R5.1 A9.0 F14.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	26
FoV read	350 mm
FoV phase	59.4 %
Slice thickness	5.00 mm
TR	27.65 ms
TE	5.51 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Elliptical filter, Image Filter
Coil elements	HEA;HEP

**Contrast - Common**

TR	27.65 ms
TE	5.51 ms
Flip angle	12 deg

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	350 mm
FoV phase	59.4 %
Slice thickness	5.00 mm
Base resolution	256
Phase resolution	65 %
Slice resolution	54 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
----------	--------

**Resolution - iPAT**

Accel. factor PE	2
Ref. lines PE	30
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	On
Intensity	Medium
Edge Enhancement	3
Smoothing	3
Unfiltered images	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R5.1 A9.0 F14.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	26
FoV read	350 mm
FoV phase	59.4 %
Slice thickness	5.00 mm
TR	27.65 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R5.1 A9.0 F14.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.1 A9.0 F14.7
R	5.1 mm
A	9.0 mm
F	14.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Special sat.	None
--------------	------

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	27.65 ms
Concatenations	1
Segments	1

**Angio - Common**

Flow mode	Free
Encodings	2
Velocity enc. 1	40 cm/s
Velocity enc. 2	40 cm/s
Direction 1	F >> H
Direction 2	Through plane
Rephased images	On
Magnitude images	Off
Magnitude sum	On
Phase images	Off

**Angio - Inline**

Subtract	Off
Measurements	1
StdDev	Off

**Angio - Inline**

Save original images	On
----------------------	----

**Angio - MIP**

MIP-Sag	On
MIP-Cor	On
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Angio - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	1
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	300 Hz/Px

**Sequence - Part 2**

Segments	1
RF pulse type	Fast
Gradient mode	Fast
RF spoiling	Off

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\PC\_Sinus

TA: 2:11 PM: FIX Voxel size: 0.7×0.7×5.0 mmPAT: 2 Rel. SNR: 1.00 : fl\_r

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	R6.3 P58.2 H77.5 mm
Orientation	C > T16.3 > S2.6
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	21.70 ms
TE	6.59 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	21.70 ms
TE	6.59 ms
TD	0 ms
Flip angle	12 deg
Wrap-up Magn.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	224
Phase resolution	100 %
Phase partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2

**Resolution - iPAT**

Ref. lines PE	32
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
POCS	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	R6.3 P58.2 H77.5 mm
Orientation	C > T16.3 > S2.6
Phase enc. dir.	R >> L
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	21.70 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R6.3 P58.2 H77.5 mm
Orientation	C > T16.3 > S2.6
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	R6.3 P58.2 H77.5
R	6.3 mm
P	58.2 mm
H	77.5 mm
Initial Rotation	0.00 deg
Initial Orientation	C > T
C > T	16.3
> S	2.6

**Geometry - Saturation**

Wrap-up Magn.	None
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R6.3 P58.2 H77.5 mm
Orientation	C > T16.3 > S2.6
Rotation	0.00 deg
R >> L	160 mm
F >> H	160 mm
A >> P	5 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	Pulse/Retro
Average cycle	No Signal ms
Average cycle	No Signal ms
Calculated phases	32
TR	21.70 ms
Concatenations	1
Segments	1
Arrhythmia detection	None

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Angio - Common**

Flow mode	Single dir.
Encodings	1
Velocity enc.	50 cm/s
Direction	Through plane
Rephased images	On
Magnitude images	On

**Angio - Common**

Magnitude sum	Off
Phase images	On

**Angio - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Angio - Cardiac**

Inline Evaluation	Off
TE	6.59 ms
TR	21.70 ms
Save original images	On

**Angio - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Angio - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Flow comp.	Yes
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	10.9 ms
Sequence type	Gre
Bandwidth	203 Hz/Px

**Sequence - Part 2**

Define	Segments
Segments	1
RF pulse type	Normal
Gradient mode	Fast*
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\t1\_fl2d\_sag\_localiser

TA: 0:51 PM: REF Voxel size: 0.7×0.7×3.0 mmPAT: Off Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	20
Dist. factor	0 %
Position	L0.9 P5.3 H70.4 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	30 %
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	150.0 ms
TE	2.46 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP

**Contrast - Common**

TR	150.0 ms
TE	2.46 ms
MTC	Off
Magn. preparation	None
Flip angle	70 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	220 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
Base resolution	320
Phase resolution	80 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	20
Dist. factor	0 %
Position	L0.9 P5.3 H70.4 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	150.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.9 P5.3 H70.4 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.9 P5.3 H70.4
L	0.9 mm
P	5.3 mm
H	70.4 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	150.0 ms
Concatenations	1
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	220 mm
FoV phase	100.0 %
Phase resolution	80 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off

**Inline - Common**

Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	70 deg
Measurements	1
Contrasts	1
TR	150.0 ms
TE	2.46 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Bandwidth	320 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	None
RF pulse type	Fast
Gradient mode	Performance
Excitation	Slice-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	60 s

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\PC\_Aqueduct

TA: 2:27 PM: FIX Voxel size: 0.6×0.6×5.0 mmPAT: 2 Rel. SNR: 1.00 : fl\_r

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	1
Dist. factor	20 %
Position	L0.6 P22.4 H72.9 mm
Orientation	T > C21.4 > S-4.8
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	22.94 ms
TE	7.06 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	22.94 ms
TE	7.06 ms
TD	0 ms
Flip angle	12 deg
Wrap-up Magn.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2

**Resolution - iPAT**

Ref. lines PE	32
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
POCS	Off

**Geometry - Common**

Slice group	1
Slices	1
Dist. factor	20 %
Position	L0.6 P22.4 H72.9 mm
Orientation	T > C21.4 > S-4.8
Phase enc. dir.	R >> L
FoV read	160 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	22.94 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.6 P22.4 H72.9 mm
Orientation	T > C21.4 > S-4.8
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L0.6 P22.4 H72.9
L	0.6 mm
P	22.4 mm
H	72.9 mm
Initial Rotation	90.20 deg
Initial Orientation	T > C
T > C	21.4
> S	-4.8

**Geometry - Saturation**

Wrap-up Magn.	None
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L0.6 P22.4 H72.9 mm
Orientation	T > C21.4 > S-4.8
Rotation	90.20 deg
R >> L	160 mm
A >> P	160 mm
F >> H	5 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Slice-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	Pulse/Retro
Average cycle	No Signal ms
Average cycle	No Signal ms
Calculated phases	32
TR	22.94 ms
Concatenations	1
Segments	1
Arrhythmia detection	None

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Angio - Common**

Flow mode	Single dir.
Encodings	1
Velocity enc.	10 cm/s
Direction	Through plane
Rephased images	On
Magnitude images	On

**Angio - Common**

Magnitude sum	Off
Phase images	On

**Angio - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Angio - Cardiac**

Inline Evaluation	Off
TE	7.06 ms
TR	22.94 ms
Save original images	On

**Angio - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Angio - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	2D
Reordering	Linear
Asymmetric echo	Weak
Flow comp.	Yes
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	11.5 ms
Sequence type	Gre
Bandwidth	201 Hz/Px

**Sequence - Part 2**

Define	Segments
Segments	1
RF pulse type	Normal
Gradient mode	Fast*
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s



## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\AAHead\_Scout\_32ch-head-coil

TA: 0:14 PM: FIX Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3.15 ms
TE	1.37 ms
Flip angle	8 deg

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Resolution - Common**

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

**Resolution - iPAT**

Reference scan mode	Integrated
---------------------	------------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off

**System - Miscellaneous**

Coil Select Mode	Default
------------------	---------

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Flip angle	8 deg
Measurements	1
Time to center	6.2 s

**Inline - Inline**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1

**Inline - MapIt**

Contrasts	1
TR	3.15 ms
TE	1.37 ms

**Sequence - Part 1**

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\IR-SPGR\_TI\_600\_v4

TA: 1:55 PM: FIX Voxel size: 1.2×1.2×1.2 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	1040.0 ms
TE	1.82 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	1040.0 ms
TE	1.82 ms
Magn. preparation	Non-sel. IR
TI	600 ms
Flip angle	5 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	1040.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.9 F2.9
L	0.0 mm
A	3.9 mm
F	2.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Rotation	97.97 deg
F >> H	192 mm
A >> P	240 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	3.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1040.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI	600 ms
Fat suppr.	None
Dark blood	Off
FoV read	192 mm
FoV phase	125.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	5 deg
Measurements	1
TR	1040.0 ms
TE	1.82 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	5.4 ms
Bandwidth	520 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On
Turbo factor	160

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\IR-SPGR\_TI\_1500\_v4

TA: 3:35 PM: FIX Voxel size: 1.2×1.2×1.2 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	1940.0 ms
TE	1.82 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	1940.0 ms
TE	1.82 ms
Magn. preparation	Non-sel. IR
TI	1500 ms
Flip angle	5 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	1940.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.9 F2.9
L	0.0 mm
A	3.9 mm
F	2.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Rotation	97.97 deg
F >> H	192 mm
A >> P	240 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	3.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	1940.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI	1500 ms
Fat suppr.	None
Dark blood	Off
FoV read	192 mm
FoV phase	125.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	5 deg
Measurements	1
TR	1940.0 ms
TE	1.82 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	5.4 ms
Bandwidth	520 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On
Turbo factor	160

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\IR-SPGR\_IROff\_2\_v4

TA: 1:36 PM: FIX Voxel size: 1.2×1.2×1.2 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	865.0 ms
TE	1.82 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	865.0 ms
TE	1.82 ms
Magn. preparation	None
Flip angle	2 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	865.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.9 F2.9
L	0.0 mm
A	3.9 mm
F	2.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Rotation	97.97 deg
F >> H	192 mm
A >> P	240 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	3.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	865.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	192 mm
FoV phase	125.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	2 deg
Measurements	1
TR	865.0 ms
TE	1.82 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	5.4 ms
Bandwidth	520 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On
Turbo factor	160

**Sequence - Assistant**

Mode	Off
------	-----



## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\IR-SPGR\_IROff\_5\_v4

TA: 1:36 PM: FIX Voxel size: 1.2×1.2×1.2 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	865.0 ms
TE	1.82 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	865.0 ms
TE	1.82 ms
Magn. preparation	None
Flip angle	5 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	865.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.9 F2.9
L	0.0 mm
A	3.9 mm
F	2.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Rotation	97.97 deg
F >> H	192 mm
A >> P	240 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	3.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	865.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	192 mm
FoV phase	125.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	5 deg
Measurements	1
TR	865.0 ms
TE	1.82 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	5.4 ms
Bandwidth	520 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On
Turbo factor	160

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b\IR-SPGR\_IROff\_12\_v4

TA: 1:36 PM: FIX Voxel size: 1.2×1.2×1.2 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	865.0 ms
TE	1.82 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	865.0 ms
TE	1.82 ms
Magn. preparation	None
Flip angle	12 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	160
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	1.20 mm
TR	865.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.9 F2.9
L	0.0 mm
A	3.9 mm
F	2.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm

**System - Miscellaneous**

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Rotation	97.97 deg
F >> H	192 mm
A >> P	240 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	3.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	865.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	192 mm
FoV phase	125.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	12 deg
Measurements	1
TR	865.0 ms
TE	1.82 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	5.4 ms
Bandwidth	520 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On
Turbo factor	160

**Sequence - Assistant**

Mode	Off
------	-----

## \\Study Protocols\BRAIN\Stroke\Investigate - E161686 v1b1t1\_vibe\_sag\_DCE\_2mm

TA: 21:08 PM: FIX Voxel size: 2.0×2.0×2.0 mmPAT: Off Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	Off
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	96
FoV read	192 mm
FoV phase	125.0 %
Slice thickness	2.0 mm
TR	3.44 ms
TE	1.68 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3.44 ms
TE	1.68 ms
Flip angle	15.0 deg
Fat suppr.	None
Water suppr.	None
Dixon	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	32
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s

**Contrast - Dynamic**

Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Multiple series	Off

**Resolution - Common**

FoV read	192 mm
FoV phase	125.0 %
Slice thickness	2.0 mm
Base resolution	96
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
POCS	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	96
FoV read	192 mm
FoV phase	125.0 %

**Geometry - Common**

Slice thickness	2.0 mm
TR	3.44 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L0.0 A3.9 F2.9
L	0.0 mm
A	3.9 mm
F	2.9 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	None
Water suppr.	None
Dixon	Off
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R8.0 A13.6 H6.2 mm
Orientation	S > T-4.2 > C-2.3
Rotation	97.97 deg
F >> H	192 mm
A >> P	240 mm
R >> L	192 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Non-sel.

**System - Tx/Rx**

Frequency 1H	123.244448 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	6.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

View sharing	Off
Flip angle	15.0 deg
Measurements	32
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Burn time-to-center	Off
Temporal interpolation	1
3D centric reordering	Off
Time to center	20.0 s

**Inline - Inline**

Subtract	Off
Measurements	32
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
---------	-----

**Inline - MIP**

MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	32
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	15.0 deg
Measurements	32
Contrasts	1
TR	3.44 ms
TE	1.68 ms

**Sequence - Part 1**

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	1
Optimization	Min. TE

**Sequence - Part 1**

Multi-slice mode	Sequential
Bandwidth	520 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On

**Sequence - Assistant**

Mode	Off
Allowed delay	0 s