

# IACAS 1

	Anatomical	Rest
Manufacturer	GE	GE
Model	Signa HDx	Signa HDx
Headcoil	8 Chan	8 Chan
Field Strength	3T	3T
Sequence	3D BRAVO	EPI
Flip Angle [Deg]	7.0	90
Inversion Time (TI) [ms]	1100	-
Echo Time (TE) [ms]	2.984	30
Repetition Time (TR) [ms]	7.788	2000
Bandwidth per Voxel (Readout) [Hz]	122	7812.5
Parallel Acquisition	ASSET x 2	-
Partial Fourier	Off	-
Number of Slices	192	32
Slice Orientation	Saggital	-
Slice Phase Encoding Direction	Right to Left	Right to Left
Slice Acquisition Order	Sequential Ascending	Interleaved Ascending
Slice Thickness [mm]	1.0	4.0
Slice Gap [mm]	0	0.6
Field of View [mm]	256	220
Acquisition Matrix	256x256	64x64
Slice In-Place Resolution [mm <sup>2</sup> ]	1.0x1.0	3.4x3.4
Number of Measurements	-	240
Acquisition Time [min:sec]	5:02	8:00
Fat Suppression	None	No
Number of Directions	-	-
Number of B Zeros	-	-
B Value(s) [s/mm <sup>2</sup> ]	-	-
Averages	-	-